

Your GLPI with Read-Only User

Medulla requires the creation of specific views in the **GLPI** database to function properly.

GLPI Read-Only

If you provide us with a user who has read-only access to your GLPI database, you will need to manually apply an SQL file before installing Medulla. This file contains the queries needed to create the required views.

Here is the link to the SQL file:

<https://dl.medulla-tech.io/nc/glpi-100.sql>

GLPI Write Access

If you allow **write** access to your GLPI database, Medulla will automatically apply the necessary views during installation.

“ These views are required for Medulla to function.

Revision #1

Created 2026-04-30 07:36:41 UTC by Adrien Thaisse

Updated 2026-04-30 07:36:41 UTC by Adrien Thaisse