

Imagination groupsg

Displays all existing imaging groups.

devdemo.siveo.lan : Page principale > Imaging > List all profiles of computers Cliquer pour passer en mode expert | ▼ Déconnexion root

- Statut du serveur d'imaging
- Gérer les masters
- Gérer les services de démarrage
- Menu de démarrage par défaut
- Scripts de post-imaging
- Liste des profils
- Ajouter un profil
- Configuration de l'imaging
- Gestionnaire de fichiers de réponse sysprep
- Tous les groupes d'imaging**
- Ajouter un groupe d'imaging

Groupes d'imaging

Elements 0 à 0 - Total 0

Nom du groupe	Favori	Actions
---------------	--------	---------

Rechercher

- Each group contains a specific configuration linked to a type of deployment (e.g. classrooms, specific computer park...)
- You can view, edit or delete them.

Add a group of imagining

Create an imagining group to organize your deployments.

medulla Dashboard kiosk Utilisateurs Groupes Ordinateurs **Imaging** Package Audit Mises à jour Sauvegarde Historique Admin

wva.medulla-tech.io : Page principale > Imaging > Computers Profile Creator Cliquer pour passer en mode expert | ▼ Déconnexion root

- Statut du serveur d'imaging
- Gérer les masters
- Gérer les services de démarrage
- Menu de démarrage par défaut
- Scripts de post-imaging
- Liste des profils
- Ajouter un profil
- Configuration de l'imaging
- Gestionnaire de fichiers de réponse sysprep
- Tous les groupes d'imaging
- Ajouter un groupe d'imaging**

Sélection du serveur d'imaging

Sélectionner un serveur d'imaging pour ce groupe

Statut du serveur d'imaging

Gérer les masters

Gérer les services de démarrage

Menu de démarrage par défaut

Scripts de post-imaging

Liste des profils

Ajouter un profil

Configuration de l'imaging

Gestionnaire de fichiers de réponse sysprep

Tous les groupes d'imaging

Ajouter un groupe d'imaging

Création d'un groupe d'imaging

Nom du groupe

Marquer en tant que favori Oui Non

Tous les ordinateurs

Membres du groupe

Steps

- Side menu **Add a group of imaging**
- Choose an entity
- Validate entity
- The group addition form appears
- Name the group
- It is possible to choose to put this group as a favorite
- Combine one or more machines
- Click **Validate**

This group will then appear in the list of available groups.

Revision #2

Created 2025-05-12 14:16:00 UTC by DocsAdmin

Updated 2025-05-28 09:59:08 UTC by DocsAdmin